Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Consultar las instrucciones de uso / Potrebo poštivati uputstva za upotrebu / Bodete ce se o upastvoto za upotrebu / Необходимо учитывать инструкции по применению / Viz návod k použití / Кодану бойынша нұсқаулықты ескеру қажет

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as the ingots have collapsed and oxide net cracks. For melting by flame rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to vibrate. Allow the cylinder to cool down slowly to the room temperature and then deflask.

**Processing:** Process the frame in one direction using a hard-metal milling tool. Remove any sharp edges.

**Ceramic:** Use ceramic, suitable for porcelain to metal restoration.

1. Sand-blast the surface by a pencil-blaster. Blast all of the surface on which the porcelain is to be applied with  $\text{Al}_2\text{O}_3$  - Interalex 250 µm and clean the surface (ultrasonically cleaning, under vapour or boiling water, diethyl acetate).

2. If oxide firing is carried out to check the surface, it need to be fired on 960 °C at vacuum for 5 - 10 minutes. Frame has to show an even dark grey colour. Oxide layer should be removed with  $\text{Al}_2\text{O}_3$  - Interalex 250 µm and cleaned again (vapour, boiling water, ethyl acetate). The opaque is then applied on the surface by a first thin wash firing followed by a second covering opaque.

3. Firing has to be done according to the ceramic manufacturer's instructions.

4. After each firing step (dentine bake, build-up and glazing) carry out slow cooling down to 700 °C.

5. Remove ceramic mechanically. Hydrofluoric acid (HF) corrodes metal frame.

**Finishing:** After firing of the metal-ceramic prosthetic work sand-blast visible surface of the frame with  $\text{Al}_2\text{O}_3$  - Interalex 50 micr. and polish it with suitable grinding and polishing instruments. At the end use Universal polishing paste for Co-Cr-Mo alloys REF 460 and polish up to high gloss.

**Soldering and Welding:** Soldering before firing of the frame can be carried out with suitable base-metal-solder and high temperature flux. For welding with laser use suitable base-metal welding wires. The requirements of EN ISO 9333 must be followed.

**Safety warnings:** Use suction unit, wear gloves, goggles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

**Reciprocal Actions!** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Side effects:** Consider allergic hypersensitivities to contents of the alloy.

**Notice:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Technical data based on our own researches and knowledge and can be used as orientation values. The user itself is responsible for appropriate use of the product.

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